



Model Number: 91965011101

Chemspec Colorfast Extraction Upholstery and Rug Cleaner, UPC 91965011101 Gallon, Color Fast

Manufacturer: Chemspec

Chemspec Colorfast Extraction Upholstery Cleaner UPC 91965011101 Gallon Color Fast

For Use In Portable Extraction Machines Designed or Adapted For Wet Cleaning Upholstery

Designed or Adapted For Wet Cleaning Upholstery

Expands the Number of Upholstery Fabrics Which Can be Wet Cleaned. Extensive research by Chemspec laboratories on the problem of reducing the number of dyes which bleed prompted the formulation of Colorfast Extraction Upholstery Cleaner. Once cleanable only through use of a dry solvent method, now many more upholstery fabrics are cleanable with the wet cleaning process using Colorfast Extraction Upholstery Cleaner. It's formulated with special dye setting agents.

An Exclusive Formulation Colorfast Extraction Upholstery Cleaner is an exclusive combination of selected cleaning agents and natural solvents formulated for safe upholstery fabric wet cleaning with extraction equipment. Colorfast Extraction Upholstery Cleaner also contains special additives which impart a soft, gentle hand to the fabric.

A Safer Solvent Colorfast Extraction Upholstery Cleaner reduces the chance of fabric dyes bleeding. Non-combustible in use dilutions, Colorfast Extraction Upholstery Cleaner provides that additional measure of safety desired during the cleaning process.

Specifically Formulated As a Cold Water Extraction Product To Prevent Color Loss When Cleaning Wool Or Any Other Fiber That Presents a Bleeding Problem
An Acid Side Detergent Formulation That Gently, Yet Effectively Cleans Even The Most Challenging Of Potential Bleeders

Controls and eliminates browning and rinses clean.

Buffered For Economy Color Fast for Rugs in concentrate has a pH of 2.87, yet when diluted maintains its acid cleaning ability with a pH of 3.2. Color Fast for Rugs allows you to clean in confidence and eliminates the need to use an acidic textile rinse. This saves you time and reduces the chance of overwetting.

Not available in California

Colorfast Extraction Upholstery Cleaner (AGHS SDS - US/CA English - 2-24-2015)

Colorfast Extraction Upholstery Cleaner 1gal (GHS Label)

Colorfast Extraction Upholstery Cleaner (Product Image - Web)

Colorfast Extraction Upholstery Cleaner (Product Image - Print)

Form:

Liquid

Equipment: Carpet Cleaning Machines > Vacuum Cleaners > HEPA Concrete Dust Slurry Hazmat Vacuums >

SKU Quantity:

Case of four 1 gal. / 3.8 L bottles

Unit Size: 1 gal. / 3.8 L bottle

Portable Extraction:

1. Mix 1 qt. (946 mL) in 1.5 gal. (5.7L) of water.
2. Clean soiled area using standard techniques for portable extraction equipment.

Truckmount Extraction with Upholstery Tool:

1. Mix 3.5 qts. (3.3L) in 5 gal. (19L) stock solution.
2. Apply evenly without over-wetting fabric on small area. Allow to dwell for 1-2 min.
3. Rinse and extract using truckmount and Upholstery Pro Tool until excess moisture is removed.
4. Repeat process until entire upholstery piece has been cleaned.

Availability: This product was added to our catalog on Sunday 06 February, 2005