



Model Number: C-C1200

Chemspec C-C1200 Liqua Gel 12/16 oz Case Canada Only Spotter

Manufacturer: Chemspec

Chemspec C-C1200 Liqua Gel 12/16 oz Case Canada Only Spotter

Liqua-Gel with Biosolv

Product Description

Effectively removes grease, tar, ink, glue, gum, adhesives, oils and other solvent-based substances. Non-toxic RTU formulation.

Available Sizes

- 16 oz bottle (LGCS - Case of 12) (US)

- 16 oz bottle (C1200 - Case of 12) (CA)

Product Specifications

- Form: Liquid concentrate

- RTU pH: 5.74 (US)

- RTU pH: 5.74 (CA)

1. Shake well before use.

2. Test in an inconspicuous area before use.

3. Apply directly to spot.

4. Agitate area with a bone spatula, then remove excess.

5. Rinse and extract.

Liqua-Gel with Biosolv (AGHS SDS - US/CA English - 7-31-2015)

Liqua-Gel with Biosolv (AGHS SDS - US/CA English - 6-03-2015)

Liqua-Gel with Biosolv (CAN) (GHS Label)

Liqua-Gel with Biosolv (GHS Label)

Liqua-Gel with Biosolv (US/EU) (Product Label)

Liqua-Gel with Biosolv (Product Image - Web)

Liqua-Gel with Biosolv (Product Image - Print)

Availability: This product was added to our catalog on Wednesday 27 January, 2016